



# CUSTOM STABLE FITTING SET MENU

Choose from these clinically proven 24 lens fitting sets (A, B or C), or customize your 12 lens set by choosing two 6 Packs out of set A, B or C.



## 14.8 Custom Stable Elite - Oblate

Designed for pediatrics, lower corneal elevations, regular corneas, small aperture/HVID

## 14.8 Custom Stable Elite - Prolate

Designed for higher cornea elevations, small aperture/HVID

## 15.8 Custom Stable Elite - Oblate

Designed for regular corneas, high astigmatism, post-surgery (refractive and corneal)

## 15.8 Custom Stable Elite - Prolate

Designed for all forms of Keratoconus and ectasias

## 16.8 Custom Stable Elite - Oblate

Designed for larger aperture/HVID, post-surgery (refractive and corneal)

## 16.8 Custom Stable Elite - Prolate

Designed for larger aperture/HVID, all forms of Keratoconus and ectasias

Set A	6 Pack of	14.8 Custom Stable Elite	Oblate
	6 Pack of	14.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Prolate
Set B	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Prolate
Set C	6 Pack of	15.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	16.8 Custom Stable Elite	Prolate
	6 Pack of	16.8 Custom Stable Elite	Oblate